

RESUSCITAIRE CHECKLIST – WRH Antenatal Ward, Meadow & TCU

Top surface	
Heater Light Clock Clean, flat working surface O2 & Air supply with pressure relief valve (replace if either < 500 psi) Suction apparatus with Mini Yankeur	DAILY CHECK TO ENSURE WORKING
Warm towels	2
Hats – large & small	1 each
Stethoscope	1
Breathing Circuit with T-Piece	1
Self-inflating BVM bag (~500ml) with blow off valve	1 (hanging on side of resus)
O2 Saturation probe & sensor	1
Top Drawer	
Masks – size 0,1	1 each
Single use, pre-assembled (in bag) Laryngoscopes – size 00, 0, 1 (sterile). (Green/metal)**DISPOSE OF AFTER USE (inc. batteries)	1 each
Endotracheal tube introducer	1
Endotracheal tubes – size 2.5, 3.0, 3.5	1 each
CO2 detector	1
Mini Yankeur/ Suction catheter – size 10	1
Neonatal Endotracheal Tube Grip	1 of each
Sorbaderm Barrier foam Applicator	1
iGel size 1 & sachet of optilube	1
Guedel Airway – size 1, 0, 00 & 000	1 of each
Bottom Drawer	
Box non-sterile gloves	1
Spare towels	2
Spare hats – large/ small	1 each
Spare breathing circuit with T-piece	1
Spare face masks size 0 and 1	1 each
Plastic bag for pre-term infant	1
Other	
Newborn life Support Algorithm	
Signature Checking List/ Spare	
NB Check list corresponds with resus number	

Emergency Drug Box Kept On Blue Paediatric Trolley